

Contents

About HyPep	1
What is HyPep?	1
License	1
Using HyPep	1
Computational requirements	1
Installation and start up	2
HyPep user interface	2
Loading data	3
Search options	6
Interpreting Results	7
Troubleshooting	7

About HyPep

What is HyPep?

HyPep is an open-source Python based algorithm designed for hybrid analysis of neuropeptide mass spectral data through sequence homology searching and accurate mass matching.

License

Using HyPep

HyPep is freely available for download from GitHub (<https://github.com/lingjunli-research/HyPep-v1.0>) and has an included user interface for increased accessibility.

Computational requirements

HyPep requires Python 3 installation (<https://www.python.org/downloads/>). The Python version for confirmed use at time of release is Python 3.10.6. HyPep is only supported for use with Windows OS, though limited functionality may be available with MacOS and Linux/UNIX. HyPep is designed for functionality with .RAW mass spectra files, a file type proprietary to Thermo Scientific instruments.

Input files are also required from companion software: PEAKS (Bioinformatics Solutions, Inc.), TopFD (Xiaowen Liu Research Group), and RawConverter (John Yates Research Group). TopFD requires .mzML input; .RAW files can be converted to .mzML file type through

MSConvert (Proteowizard) and can be converted to .MS2 file type through RawConverter or MSConvert.

HyPep can be run using the Windows command prompt, but it is generally a good practice to run through an Anaconda environment.


Software package	Version verified for use	Download location
Python	3.10.6	https://www.python.org/downloads/
PEAKS	PEAKS Studio Xpro	https://www.bioinform.com/peaks-studio/
RawConverter	1.2.0.1	http://fields.scripps.edu/rawconv/
TopFD	1.5	https://www.toppic.org/software/toppic/register.html
MSConvert	3	https://proteowizard.sourceforge.io/download.html
Anaconda	3.9	https://www.anaconda.com/products/distribution

Installation and start up

1. Download HyPep from GitHub (<https://github.com/lingjunli-research/HyPep-v1.0>), save to any location in the C: drive.
2. Open command prompt (either Anaconda or built-in)
3. Navigate to downloaded HyPep folder within the command prompt
4. Launch HyPep GUI using the input: `python gui_flattened.py`


HyPep user interface

The HyPep user interface has three primary pages: Add Modules, Select Files, and Run Analysis, designed to access in this order.



Each parameter of Add Modules is described in later sections. It should be noted that prior to moving to the Select Files page, “Save Preferences” must be selected.

Each selection of Select Files is described in later sections. It should be noted that prior to moving to the Run Analysis page, “Save Preferences” must be selected



After the Add Modules and Select Files pages have been completed, “Run Analysis” within the Run Analysis page can be selected. Analysis has begun when the green bar begins to move. All progress of analysis will be reported in the command prompt.

Loading data

In selecting data, simply click the “browse” button to the right of the input field within the Select Files page. Unless noted, all .csv/.txt file headers must match the template files.

Sample name

Input a sample name within this field, no spaces or special characters are permitted.

The database is input as a .csv file with 3 columns: Accession, Sequence, Monoisotopic $[M+H]^+$ m/z value.

	A	B	C	D
1	Accession	Sequence	Monoisotopic Mass	
2	1	LRVAPEESPVL	1209.68376	
3	2	LRVAPEEHPVLL	1372.7947	
4	3	YAFGL(Amidated)	569.30821	
5	4	YSFGL(Amidated)	585.30312	
6	5	YTFGL(Amidated)	599.3188	
7	6	YDFGL(Amidated)	613.298	
8	7	SYAFGL(Amidated)	656.3402	
9	8	SYDFGL(Amidated)	789.3361	

The HyPep input file from PEAKS is the *de novo peptides.csv* file.

A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	
1	Fraction	Source File	Feature	Peptide	Scan	Tag Length	Denovo Score	ALC (%)	length	m/z	z	RT	Predict RT	Area	Mass	ppm	PTM	local confidence (%)	tag ≥0%	mode
2	27.2021	0817_Brain_1.raw	F2748539	ESDAAEV	F2716840	14	99	99	14	753.717	2	34.24	1.17E+07	1504.726	2.2			100 100 100 100 100	1	ESDAAEVHK HCD
3	27.2021	0817_Brain_1.raw	F2748996	ESDAAEV	F2725238	27	99	99	27	759.1259	4	52.92	3.86E+07	3032.467	2.4			100 100 100 100 100	1	ESDAAEVHK HCD
4	27.2021	0817_Brain_1.raw	F2743669	ESDAAEV	F2719200	19	99	99	19	697.6764	3	39.64	1.46E+07	2090.002	2.8			100 100 100 100 100	1	ESDAAEVHK HCD
5	27.2021	0817_Brain_1.raw	F2716596	EVGKLLK	F2728002	8	99	99	8	450.3006	2	58.71	1.92E+06	888.5851	1.7			100 100 100 100 100	1	EVGKLLK HCD
6	27.2021	0817_Brain_1.raw	F2720935	KSVFDQLK	F2721488	8	99	99	8	482.7776	2	44.63	1.42E+06	963.5389	1.9			100 100 100 100 100	1	KSVFDQLK HCD
7	27.2021	0817_Brain_1.raw	F2728402	EAEATLK	F2729260	7	99	99	7	381.2062	2	17.68	4.27E+06	760.3967	1.5			100 100 100 100 100	1	EAEATLK HCD
8	27.2021	0817_Brain_1.raw	F2728402	EAEATLK	F2729266	7	99	99	7	381.2062	2	17.68	4.27E+06	760.3967	1.5			100 100 100 100 100	1	EAEATLK HCD
9	27.2021	0817_Brain_1.raw	F2735872	TMDVYK	F2712322	10	99	99	10	611.7941	2	24.15	1.72E+06	1221.57	2.9			99 100 100 100 100	10	TMDVYK HCD
10	27.2021	0817_Brain_1.raw	F2720935	KSVFDQLK	F2721493	8	99	99	8	482.7776	2	44.63	1.42E+06	963.5389	1.9			100 100 100 100 100	1	KSVFDQLK HCD
11	27.2021	0817_Brain_1.raw	F2716983	KFDPALK	F2713180	8	99	99	8	453.2583	2	25.93	7.02E+06	904.5018	0.3			100 100 100 100 99	10	KFDPALK HCD
12	27.2021	0817_Brain_1.raw	F2710230	RVEVSYSR	F2726922	9	99	99	9	398.5516	3	11.75	3.39E+06	1000.100	1.6			99 99 100 100 100	10	RVEVSYSR HCD
13	27.2021	0817_Brain_1.raw	F2748539	ESDAAEV	F2716852	14	99	99	14	753.717	2	34.24	1.17E+07	1504.726	2.2			99 98 100 100 100	10	ESDAAEVHK HCD
14	27.2021	0817_Brain_1.raw	F2716983	KFDPALK	F2713122	8	99	99	8	453.2583	2	25.93	7.02E+06	904.5018	0.3			100 100 100 100 99	10	KFDPALK HCD
15	27.2021	0817_Brain_1.raw	F279899	LGPMMLK	F2733200	7	99	99	7	395.2129	2	69.73	3.34E+06	788.4288	0.5			99 100 100 100 100	10	LGPMMLK HCD
16	27.2021	0817_Brain_1.raw	F2735409	EAFATIK	F2730934	11	99	99	11	607.8369	2	65.09	2.27E+08	121.655	3.3			99 99 100 100 100	10	EAFATIK HCD

The HyPep input file from TopFD file is a .csv file with a title ending in “frac.mzrt” after output from TopFD.

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
1	ID	Fraction	Envelope	Mass	MonoMz	Charge	Intensity	mzLo	mzHi	rtLo	rtHi	color	opacity	promex_score	
2	0	0	14	418.259	419.266	1	1.75E+10	419.166	421.371	99.2908	99.4825	##FF0000	0.1	-339.146	
3	1	0	6	408.276	409.284	1	6.94E+09	409.184	411.389	99.0657	99.1578	##FF0000	0.1	-339.146	
4	2	0	3	5608.91	1403.23	4	1.52E+08	1403.13	1405.59	98.488	98.5062	##FF0000	0.1	53.0046	
5	2	0	11	5608.91	1122.79	5	2.85E+09	1122.69	1124.69	98.4687	98.5265	##FF0000	0.1	53.0046	
6	2	0	16	5608.91	935.826	6	6.04E+09	935.726	937.429	98.4298	98.5448	##FF0000	0.1	53.0046	
7	2	0	16	5608.91	802.28	7	1.08E+10	802.18	803.669	98.4388	98.5448	##FF0000	0.1	53.0046	
8	2	0	12	5608.91	702.121	8	3.17E+09	702.021	703.349	98.4298	98.5665	##FF0000	0.1	53.0046	
9	3	0	11	1182.77	592.392	2	7.59E+08	592.292	593.996	9.91483	10.0148	##FF0000	0.1	54.4977	
10	3	0	18	1182.77	395.264	3	9.58E+09	395.164	396.367	9.893	10.132	##FF0000	0.1	54.4977	
11	3	0	11	1182.77	296.7	4	2.78E+09	296.6	297.552	9.91483	10.0488	##FF0000	0.1	54.4977	
12	3	0	1	1182.77	237.561	5	4.74E+06	237.461	238.263	9.927	9.927	##FF0000	0.1	54.4977	

The HyPep input file from RawConverter or MSConvert is a .MS2 file.

```


2021_0817_Brain_1 - Notepad
File Edit Format View Help
H Creation Date 8/1/2022 9:43:21 AM
H Extractor RawConverter
H ExtractorVersion 1.1.0.18
H Comments RawConverter written by Lin He, 2014
H Comments RawConverter modified by Yen-Yin Chu, 2015
H Comments RawConverter modified by Rohan Rampuria, 2016
H ExtractorOptions MSn
H AcquisitionMethod Data-Dependent
H InstrumentType FTMS
H DataType Centroid
H ScanType MS2
H Resolution
H IsolationWindow
H FirstScan 1
H LastScan 56147
H MonoIsotopic PrecMz False
S 000002 000002 350.95926
I RetTime 0.01
I IonInjectionTime 119.999997317791
I ActivationType HCD
I InstrumentType FTMS
I TemperatureFTAnalyzer -1
I Filter FTMS + c NSI d Full ms2 350.9593@hcd30.00 [100.0000-735.1
I PrecursorScan 1
I PrecursorInt 64845.1
Z 2 700.91124
101.0599 1361.1 0 15800
104.1071 761.9 0 14400
105.8739 645.1 0 13900
106.0521 772.1 0 14900
114.3627 617.3 0 15100
116.986 1298 0 15500
123.0552 2568.6 0 20800
126.055 2993.8 0 19100
136.0617 2949.8 1 19004
137.0656 644.5 1 13700
141.8851 719.6 0 12700
146.311 837.4 0 14200
147.7578 811.8 0 13500
153.9545 628 0 12200
154.0497 3277 0 17900
156.9298 641.6 0 11100
Ln 1, Col 1 100% Unix (LF) UTF-8

```

ProteinProspector Input

A comma-delimited .txt file is required for each database sequence, which contains the ion mass, ion identity, and the ion chemical formula. If ion chemical formula is not known, keep the column header and use this column to repeat the ion identity. Column headers should be: mass, ion, formula. The file name for each sequence should begin with "Theoretical_b_y_ion_" followed by the peptide name, written exactly as it is written in the database. These ion lists can be generated from ProteinProspector (<https://prospector.ucsf.edu/prospector/cgi-bin/msform.cgi?form=msproduct>), or another program of choice.

*An alpha version of HyPep 2.0 is available upon request, which has this function automated.



```
mass,ion,formula
70.0651,P,C4 H8 N1
86.0964,L,C5 H12 N1
102.055,E,C4 H8 N1 O2
120.0808,F,C8 H10 N1
126.055,P,C6 H8 N1 O2
131.1179,y1,C6 H15 N2 O1
136.0757,Y,C8 H10 N1 O1
182.0812,a2,C9 H12 N1 O3
188.1394,y2,C8 H18 N3 O2
192.0655,b2-H20,C10 H10 N1 O3
210.0761,b2,C10 H12 N1 O4
335.2078,y3,C17 H27 N4 O3
345.1445,a3,C18 H21 N2 O5
355.1288,b3-H20,C19 H19 N2 O5
373.1394,b3,C19 H21 N2 O6
406.2449,y4,C20 H32 N5 O4
416.1816,a4,C21 H26 N3 O6
426.166,b4-H20,C22 H24 N3 O6
444.1765,b4,C22 H26 N3 O7
563.25,a5,C30 H35 N4 O7
569.3082,y5,C29 H41 N6 O6
573.2344,b5-H20,C31 H33 N4 O7
591.2449,b5,C31 H35 N4 O8
620.2715,a6,C32 H38 N5 O8
630.2558,b6-H20,C33 H36 N5 O8
648.2664,b6,C33 H38 N5 O9
666.361,y6,C34 H48 N7 O7
760.3665,MH-H20,C39 H50 N7 O9
778.377,MH,C39 H52 N7 O10
```

Output Directory

The output directory can be located in any drive and is an empty folder where all results will be exported. If folder is not empty, previous results will be written over. The output folder must not be a sub-directory of the input folder.

Search options

False discovery rate algorithm selection

Reverse

Reverses each database sequence in original target database.

Shuffle

Randomly shuffles amino acid residues within each sequence in original target database.

Random

All amino acids in original target database are concatenated into one sequence, shuffled, and new sequences are formed while maintaining the sequence length frequency from the original target database.

Hybrid

A novel decoy database generation method which is a hybrid of shuffle and random decoy database methods. All amino acids in original target database are concatenated into one sequence, shuffled, and new sequences are formed while 1) maintaining the sequence length frequency from the original target database and 2) the amino acid frequency pattern of each sequence from the original target database is conserved in each corresponding decoy sequence.

False discovery rate threshold

The false discovery rate corresponds to the expected number of theoretical false positive identifications.

Sliding window size

The sliding window size controls the size of the amino acid window during alignment between a *de novo* sequenced peptide query and database sequence. The sliding window size ranges from 1-10. Must be a whole number. Recommended value is 2.

Maximum precursor charge

The maximum precursor charge is the maximum charge value that will be considered for a precursor peak. All charges from 1- n , where n is the maximum precursor charge, will be analyzed. Must be whole number. Recommended value is 8.

Maximum fragment charge

The maximum fragment charge is the maximum charge value that will be considered for a fragment peak. All charges from 1- n , where n is the maximum fragment charge, will be analyzed. Must be whole number. Recommended value is 4.

Precursor error cutoff

Error margin for a precursor match.

Fragment error cutoff

Error margin for a fragment match.

Promex score cutoff

Promex is a likeliness score exported from TopFD, where higher values represent greater likeliness. Recommended value is -10.

Number of matching loops

The number of loops that will be conducted to make PSMs while maximizing the number of unique peptide identifications (*i.e.*, more loops will result in more identifications to an extent). Recommended value is 5.

Interpreting Results

All results will be located in the specified output folder, in the sub-folder named “final_results”. Results from the SHS module are called “SampleName_SHS_final_report”, the results from the AMM module are called “SampleName_AMM_final_report”, and the final HyPep output containing the combined SHS and AMM module results is labeled “SampleName_combined_AMM_SHS_final_report”. Additionally, all parameters selected prior to analysis are reported as a parameter .txt file.

Troubleshooting

With any issues, please either submit a public issue on the HyPep GitHub (<https://github.com/lingjunli-research/HyPep-v1.0/issues>), or email lawashburn@wisc.edu with questions. Future releases of HyPep will be announced at <https://www.lilabs.org/resources>.